

*Effective technologies designed
for decontamination.*



Hydrogen Peroxide solution designed to
disinfect efficiently and very quickly all
the places where hygiene is essential

Medtra commercializes and
develops products as well as
technology which act
effectively on any kind of
infection with the aim of
guaranteeing the best
possible hygiene for
everyone.

**Medtra technical,
H₂O₂ solution
concentrated at 7%,
provides effective
and fast disinfection,
spray safely**

Medical Device : 0499

ISO 13485 /9001

Technical informations

Usage	Ready to use solution
Actif ingrédient	Hydrogen peroxyde 7%
GMDN Nr.	45058
Parfumes	2 flavors available: Neutral, apple
Water solubility	Total
Shelf life	17 months unopened 2 months after opening
Chemical stability	Stable from 5°C to 30°C
Relative density	1
Storage conditions	In original packaging Protect from sunlight Avoid temperature fluctuation
Packaging	Bottles of 1L (boxes of 6 bottles) Can of 5L

Efficiency

Tests were carried out on Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, Aspergillus Niger, Penicillium SPP, Norovirus, Bacillus stearothermophilus, ...

Applications

The Medtrapro technical is applicable in various fields:

hospitals, nursing homes, ambulances, dental offices, pharmaceutical, para-medical community

Warnings

- During any manipulation, wear an appropriate protection outfit (gloves, eye protection, ...)
- Check the expiry date of the product
- Read the safety sheet (available on our website: www.medtraworld.com)
- Handle the product away from any source of heat

Disinfection by air is an additional tool for mechanical cleaning methods and must be in accordance with the cleaning protocol applied in the institution

Benefits

- ✓ Quick disinfection of the volumes and areas in less than 1 hour
- ✓ Ready to use solution
- ✓ No residue
- ✓ No residual humidity
- ✓ The product (H₂O₂) is ecologically inoffensive and biodegradable
- No corrosion on most of the areas

